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




Use instructions

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KaVo UniK

Fabricante:

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CONGRATULATIONS!

Dear Doctor.

You have just acquired the **KaVo Unik**.

A certain satisfaction and reliability only a world wide recognized company like KaVo is able to offer to you.

For you to work safely and economically, it is essential to observe all recommendations in this manual. This will provide a good performance, reliability and benefits from your device.

To reproduce or provide third parties with these instructions, prior permission from KaVo do Brasil Indústria e Comércio Ltda is required.

The KaVo Group is always researching new technologies and innovations to improve their products. The technical specifications of this product and the User's Manual might change based on new improvements, without a prior notice.

We inform all possible modifications does not results in right to retool any existing device.

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ÍNDIX

USE	6
1 APPLICATION	6
2 PRINCIPLE OF FUNCTIONING	6
3 IMPORTANT NOTES	6
4 IMPROPER USE	12
5 PRODUCT PRESENTATION	14
6 DESCRIPTION	14
7 OPERATION	18
7.1 Main switch	18
7.2 Adjustment of headrest height	18
7.3 Adjustment of headrest inclination	18
7.4 Equipment table column height adjustment	18
7.5 Adjustment of table height	18
7.6 Spray water reservoir	18
7.7 Chair operation	20
7.8 Operating the dentist's element	22
7.8.1 Coupling the handpiece on the tubing	22
7.8.2 Spray preselection	22
7.8.3 Spray regulation	22
7.8.4 Selectiong the optic fiber light	22
7.8.5 Use the X-ray-viewer	22
7.8.6 Handpiece driving	24
7.8.7 Handling of handpieces	24
7.8.8 3F handpiece	24
7.8.9 Using the SONICborden	24
7.8.10 Using the ULTRA scaler	24
7.9 Operating the suction unit / assistant's element	24
7.9.1 Handling of suction tubings	24
7.9.2 Using the Poly 600 (accessory)	26
7.9.3 Water jet for spittoon cleaning, through the register	26
7.9.4 Water jet for spittoon cleaning, with timed interruption, activated in the foot control Joystick	26
7.9.5 Supply of the glass (accessory)	26
7.9.6 Movement of the suction unit	26
7.10 Operating the operating light	26
8 CLEANING / DISINFECTION	30
8.1 Table with guidelines for routine preservation	30
8.2 Technical safety controls	32
8.3 Annualy inspections	34
8.4 Assepto Sys system (accessory)	34
9 FAILURE ELIMINATION	36
10 OPTIONS AND ADAPTABLE ACCESSORIES	38
11 SIMBOLOGY	44
12 TECHNICAL CHARACTERISTICS	46
13 INSTALLATION	48
14 WARRANTY TERM	50

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USE

1 APPLICATION

KaVo UniK comprises a chair designed to accommodate the element "patient", an equipment to drive and control 3-in-1 syringes, micro motors, turbines and other instruments used in the dental practice, a suction unit to collect debris by means of the spittoon/ejectors and an operating light to lighten the oral cavity, complying with NBR ISO 6875 and NBR IEC 60 601-1 norms.

It may only be used by a competent person for the described application, observing the valid accident prevention measures in the present operating instructions.

It is the duty of the user to use only faultless tools, to ensure that they are used for the correct purpose and to protect patients and third parties from all dangers.

2 PRINCIPLE OF FUNCTIONING

Unit supplied by voltages of 110, 127 or 220 – 50/60HZ, for which it requires adaptation in the output connections of the transformer. Motor reducers with electro-electronic activation (using a micro-controlled system), supplied with 24 VDC and controlled via the foot control and membrane – micro-interruptors), are responsible for moving the backrest and the seat of the chair.

The suction unit is fixed to the chair and has a rinsing bowl, connected to the draining network and a venturi type suction system (injector activated by compressed air) for 6.3mm nozzles. A second venturi suction system for 9.5 mm nozzles and a 3-function syringe (air, water and spray) can be incorporated as accessories. The pressurized water tanks for the spray of the handpieces and the Aspetto Sys (asepsis system of the internal channels of the equipment handpiece hoses, except for the triple syringe) are located in the suction system.

The equipment has control blocks for two handpieces, fed with air supplied by the network and water provided by the pressurized tank. Controlled by an outflow regulating valve located in the foot control, they execute the activation and rotation control of the handpieces. A control block for a third handpiece can be incorporated as an accessory.

Appliance of oral illumination of intensity single, supplied with 24 VAC coming from the transformer and controlled by a micro-interruptor located in the foot control.

3 IMPORTANT NOTES

KaVo will not assume liability for damage caused by:

- External effects (poor quality of the power supply or incorrect installation of this).
- Use of information provided by third parties.
- Repairs made by personnel unauthorized by KaVo.
- Use of accessories, spare parts and materials not specified by KaVo do Brasil.

KaVo products may be repaired and serviced only by technicians specially trained and authorized by KaVo.

Only original spare parts may be used to operate and repair the equipment.

⚠ KaVo of Brazil does not take any responsibility for damages caused by use of accessories, spare parts and materials not specified by it.

Statutory provisions

The relevant general directives/guidelines and/or national law, national regulations and technical rules are to be applied and fulfilled for starting up and operating the KaVo product for the intended purpose.

Biocompatibility

Parts of this equipment and its accessories which get in contact with the human body are according to the directives and principles of the norm ISO 10993-1.

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Normative classification of the product

- According to the kind of protection against electrical shocks: Class I Equipment
- According to the grade of protection against electrical shocks: Odontological chair of the patient and Operating light: applied part type B
Dental unit: applied part type BF
- According to the grade of protection against dangerous water penetration: IPX0
- According to the grade of application safety in the presence of an inflammable anesthetic mixture: Equipment not suitable for use in the presence of an inflammable anesthetic mixture with air, oxygen or nitrous oxide
- According to the function mode:
Odontological chair of the patient: continuous function with intermittent load. The time of operation is of 20 seconds and the time of pause is of 200 seconds.
Operating light and dental unit: continuous function.

Electromagnetic compatibility

Instructions an declaration of the manufacturer – Electromagnetic emissions

The KaVo Unik Unit requires special precautions as regards its Electromagnetic Compatibility and needs to be installed and placed in functioning as per the information about Electromagnetic Compatibility provided in this Manual. In general, the mobile and portable RF communication equipment can affect electromedical equipment.

⚠ The use of accessories and options other than those specified by KaVo may result in an increase of emissions or decrease of the immunity of the equipment or system. All the accessories and options with which the equipment was tested are described in this manual. This equipment is intended to be used in an electromagnetic environment as specified below.

The user of this equipment should ensure that it is used in this type of environment.

Emission Tests	Compliance	Electromagnetic Environment – directives
RF emissions CISPR 11	Group 1	The KaVo Unik Unit does not only use RF energy for its internal functions. Nevertheless, its emissions of RF are very low and it is unlikely that they cause any interference in electronic equipment which is nearby.
RF emissions CISPR 11	Class A	The KaVo Unik Unit suitable to be used in all establishments, including residential establishments and those connected directly to the public distribution network of low voltage electricity which supplies buildings for domestic use, provided that the installation recommendations contained in the installation manual are obeyed.
Emissions of Harmonics IEC 61000-3-2	Not applicable	
Emissions owing to the fluctuation of voltage / flickering IEC 61000-3-3	Not applicable	

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Directives and declaration of the manufacturer - electromagnetic immunity

The KaVo UniK Unit is intended to be used in the electromagnetic environment specified below.

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Immunity test	Test Level NBR IEC 60601	Level of Compliance	Electromagnetic Environment - Directive
RF Conducted IEC 61000-4-6	3 Vrms 150 kHz up to 80 MHz	[3] V	Portable and mobile RF communication equipment should not be used near any part of the KaVo UniK Unit, including cables, with a separation distance less than that recommended, obtained by means of the equation applicable to the frequency of the transmitter. Recommended Distance of Separation $d = 2\sqrt{P}$
RF Radiated IEC 61000-4-3	3 V/m 80 MHz up to 2,5 GHz	[3] V/m	$d = 2\sqrt{P}$ 80 MHz up to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz up to 2.5 GHz Where P is the maximum nominal output power in watts (W), according to the manufacturer of the transmitter, and d is the recommended distance of separation in meters (m). The field intensity established by the RF transmitter, as determined by means of an electromagnetic inspection in location a, should be less than the compliance level in each frequency range b. Interference may occur around the equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the greatest frequency range is applied.

Note 2: These directives may not be applicable in all the situations. The electromagnetic propagation is affected by the absorption and reflecting of structures, objects and people.

^a The field intensities established by the fixed transmitters, such as base radio stations, telephones (cell phone/wireless), regions of mobile radios, amateur radio, AM and FM radio transmission and TV transmission cannot be foreseen theoretically with accuracy. In order to evaluate the electromagnetic environment due to fixed RF transmitters, an electromagnetic inspection of the region must be considered. If the average of the field intensity in the location in which the KaVo UniK unit is used exceeds the compliance level employed above, the KaVo UniK unit should be inspected in order to verify if it is operating normally. If any abnormality is noted, additional measures may be necessary, such as reorienting or repositioning the KaVo UniK unit.

^b Above the frequency range 150 kHz through 80 MHz, the field intensity should be less than 3 V/m.

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Immunity Tests	Test Level NBR IEC 60601	Level of Compliance	Electromagnetic Environment – directives
Electrostatic discharge (DES) as per IEC 61000 -4-2	Direct contact ± 6 kV Discharge by the air ± 8 kV	Direct contact ± 6 kV Discharge by the air ± 8 kV	The floor used for installing this equipment should be made of wood or concrete or ceramic tiles. If the floor is made of synthetic material, the relative humidity in the environment should be at least 30%.
Quick electric transitions / bursts as per IEC 61000 -4-4	Supply line ± 2 kV Communication line ± 1 kV	2	The quality of the power supply of this equipment should meet the prescription of norms.
Bursts in the supply line as per IEC 61000 -4-5	Phase - Phase ± 1 kV Phase - GND ± 2 kV	3	The quality of the power supply of this equipment should meet the prescription of norms.
Drops of voltage and brief interruptions and variances of voltage IEC 61000 -4-11	<5% UT 40% UT 70% UT <5% UT	<5% UT 40% UT 70% UT <5% UT	In general, brief interruptions and variances of voltage can affect electromedical equipment. If there are problems in the network, turn the equipment off and wait for the normal re-establishment of the power.
Magnetic Fields (50/60Hz) as per IEC 61000 -4-8	3 A/m	3 A/m	It is recommended that this equipment is not installed very close to intense magnetic fields. If this is necessary, it is advisable to verify its full operation before using it.

Environmental protection

Owing to the extensive range of products generating residue which there are in the market, their use in dental practice, the new discoveries which occur every day and the need to prevent and reduce the risks to health and the environment by correctly handling this residue, KaVo asks the user to see the health services, legal requirements/norms and directions of the manufacturers of these products to verify the appropriate procedures for determining the end destination of each type of residue generated by the products used together with this equipment.

4 IMPROPER USE

The maximum load of 135 kg for the upward movement must not be exceeded.
Load distribution (NBR ISO 6875)

Parts	Mass distribution (kg)
Head and neck	10
Superior part of the body and arms	45
Lower part of the body and arms. Hands. Superior part of the legs.	55
Legs and feet	25
TOTAL	135

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KaVo UniK

Support materials

- Use Instructions
- Authorized Technical Service (relation)
- Delivery Protocol

6 DESCRIPTION

Chair

- 1 Headrest
- 2 Backrest
- 3 Foot control
- 4 Chair base
- 5 Main switch
- 6 Fuse (inlet net)
- 7 Seat

Foot control

- 8 Key intended for activating/controlling the speed of the handpieces
- 9 On/off key and light intensity regulation of the mouth lighting apparatus
- 10 Key intended for selecting automatic functions and programming of the working positions.
- 11 Lever for activating the timing system of water for the rinsing bowl.
- 12 Standard - Activating buttons the movements of the seat/backrest.
Joystick - Plate activating the movements of the seat / backrest, the automatic movements when selected by means of key 10 and interruption of the automatic movement.

Suction unit / assistant's element

- 13 Spittoon bowl water flow control register
- 14 Water cock for the glass (accessory)
- 15 Driver for glass water (accessory)
- 16 Water spout for spittoon bowl cleaning
- 17 Spittoon bowl
- 18 Saliva Suction tubing 6,3 mm
- 19 Suction tubing 9,5mm (accessory)
- 20 Assepto Sys reservoir (accessory)
- 21 Spray water reservoir
- 22 Suction system filters (accessory)
- 23 3F handpiece holder (accessory)
- 24 Poly 600 – photopolymerizer of resins (accessory)
- 25 Time regulator of the cleaning jet in the rinsing bowl, in the timed system

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



Dentist's element

- 26 Objects tray (accessory)
- 27 Handle
- 28 3F handpiece
- 29 Handpiece holder
- 30 Screw for regulation flow (spray air)
- 31 Screw for regulation flow (spray water)
- 32 Assepto Sys system driver (accessory)
- 33 Optic fiber on/off switch (accessory)
- 34 X-ray viewer on/off switch (accessory)
- 35 X-ray viewer (accessory)
- 36 Pneumatic brake system command (optional)
- 37 ULTRAscaler frequency selector (accessory)
- 38 Viewer of the luminous intensity of the operating light.

Operating light

- 39 Handle
- 40 Mirror
- 41 Mirror protection
- 42 Lamp
- 43 Lamp protection

Dentist's element keyboard

- 44 Optic fiber keyboard (accessory).
- 45 Key for water jet commander to supply water in the cup. The action of this key in the touchpad only will be available when water set is in the cup.
- 46 ON-OFF key of the operating light device and for adjustment of light intensity, and choosing from of the selection mode
- 47 Key for activation of the movement back to the last automatic position
- 48 Light intensity Indicator
- 49 Key for activation of the chair upward movement
- 50 Key for activation of the backrest forward movement
- 51 Key for activation of the backrest rearward movement
- 52 Key for activation of the chair downward movement
- 53 Key for activation of the movement of the in-out of patient position and for setting of this position
- 54 Key for activation of the movement to the work position 1 and setting of this position
- 55 Key for activation of the movement to the work position 2 and setting of this position
- 56 Key to action of water in the bowl. The action of this key only will be available when Temporizer set of the water is in the bowl.

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
7 OPERATION

7.1 Main switch

The general switch (5) is located at the distribution box and it energizes/deenergizes the whole set.

7.2 Adjustment of headrest height

Adjustment of headrest height is made by sliding the rods through the headrest guide.

7.3 Adjustment of headrest inclination

Standard headrest

Headrest inclinations is done by actuating button (62) and turning it in the direction of the arrow.

Double-jointed headrest (accessory)

Make the adjustment by turning the button (62) counter-clockwise, placing the headrest in the desired position and turning the button (62) to the right.

7.4 Equipment table column height adjustment

Raise the coupling's protection. Take the rotating arm out. Open the plastic bushing (63) and position it on the channel corresponding to the desired height.

7.5 Adjustment of table height

Activate the command (36) of the brake, under the equipment table. Position the table at the desired height. Deactivate the command.

Movement with the brake on damages the table braking system.

Don't pull the handpieces tubings to move the table.

7.6 Spray water reservoir

Check the level of the spray tank and fill up if necessary.

To fill up:

- unscrew the reservoir (21);
- fill it up with filtered water;
- thread the reservoir (bottle) to its connection again only with the hands, enough to lean the bottleneck on the sealing rubber. After putting back the reservoir, the heard noise is not caused by air leakage, but by the pressurisation of the reservoir.

The excessive tightening causes damage to the sealing set.

After refilling and putting back the reservoir, drive an handpiece with spray for some seconds to allow the air to go out the tubings.

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
7.7 Chair operation

That happens to the movement by foot control, all instruments should be placed in brackets (table for 5T).

Chair up

Standard - Press button (12 - A);
Joystick - Actuate plate (12) in the direction **A**;
Keyboard - Press (49).

Chair down

Standard - Press button (12 - B);
Joystick - Actuate plate (12) in the direction **B**;
Keyboard - Press (52).

Backrest forward

Standard - Press button (12 - C);
Joystick - Actuate plate (12) in the direction **C**;
Keyboard - Press (50).

Backrest rearward

Standard - Press button (12 - D);
Joystick - Actuate plate (12) in the direction **D**;
Keyboard - Press (51).

Position in-out-of patient

- Standard - Briefly press button (10) and press button 12 - B.
- Joystick - Briefly press button (10) and actuate plate (12) in the direction **B**.

- In Keyboard press the key (53).

Both seat and backrest will be automatically moved to a position in-out-of patient.

This position may be changed as follow:

- Manually move the chair to the desired height and the backrest to the desired tilt.

- Press button (10) or (53) during 5-6 seconds.

- Press button 12 - B (Standard); Actuate plate (12) in the direction **B**. To confirm the recording, a beep will sound and the lamp in the lighting unit will blink twice.

The new parameters will be stored in the equipment memory.

Automatic movement to a first preset position

- Standard - Briefly press button (10) and press button 12 - A.
- Joystick - Briefly press button (10) and actuate plate (12) in the direction **A**.

- In Keyboard press the key (54).

Both seat and backrest will be automatically moved to a position in-out-of patient.

This position may be changed as follow:

- Manually move the chair to the desired height and the backrest to the desired tilt.

- Press button (10) or (54) during 5-6 seconds.

- Press button 12 - A (Standard); Actuate plate (12) in the direction **A**. To confirm the recording, a beep will sound and the lamp in the lighting unit will blink twice.

The new parameters will be stored in the equipment memory.

Automatic movement to a second preset position

- Standard - Briefly press button (10) and press button 12 - D.
- Joystick - Briefly press button (10) and actuate plate (12) in the direction **D**.

- In Keyboard press the key (55).

Both seat and backrest will be automatically moved to a position in-out-of patient.

This position may be changed as follow:

- Manually move the chair to the desired height and the backrest to the desired tilt.

- Press button (10) or (55) during 5-6 seconds.

- Press button 12 - D (Standard); Actuate plate (12) in the direction **D**. To confirm the recording, a beep will sound and the lamp in the lighting unit will blink twice.

The new parameters will be stored in the equipment memory.

Automatic movement to the last position

- Standard - Briefly press button (10) and press button 12 - C.
- Joystick - Briefly press button (10) and actuate plate (12) in the direction **C**.

LP

- In Keyboard press the key (47).

The last movement to be undone.

Interruption of the automatic movement

To interrupt the automatic movement press any of the keys from Keyboard, press any button (12) or actuate plate (12) in any direction.

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7.8 Operating the dentist's element

7.8.1 Coupling the handpiece on the tubing

On Borden and Midwest tubings, check the coincidence of the handpieces tubes with the holes on the tubing, shrink them on and thread the fixing nut only with the fingers.

On Multiflex LUX tubings, shrink the handpiece on the Multiflex coupling and tighten it against the coupling until hearing the "click".

7.8.2 Spray preselection

To select the desired status (with or without spray), open or shut off the faucet on the hose. The mark alignment indicates an open faucet.

7.8.3 Spray regulation

The spray's air/water volume can be regulated by means of the screws respectively, under the equipment's table.

(30): air
(31): water

The screws should not be completely unscrewed. The others screws existing at that same place should be adjusted by an authorized technician only. A misplaced adjustment may cause the handpiece wearing due to excessive pressure.

7.8.4 Selectiong the optic fiber light

To select the desired state (with light or without light), press the on-off switch (33) and (44) located under the equipment's table.

7.8.5 Using the X-ray-viewer

To switch the x-ray-viewer on/off, press key (34).

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
Para ligar-desligar o negatoscópio acionar a tecla (34).

7.8.6 Handpiece driving

Remove the handpiece from its support, and press the key (8) on the foot control. Rotation will vary according to pressure exerted on the key.

7.8.7 Handling of handpieces

To handle handpieces (turbines and air motors) observe operating instructions that accompany them.

7.8.8 3F handpiece

- For safety reasons only use KaVo original nozzles.
- Fit nozzle and check for correct coupling.
- Check air and water passage through the nozzle before using on the patient.

Handling

Remove 3F handpiece from holder.

Air: press key (64)

Water: press key (65)

Spray: press keys (64) and (65)

7.8.9 Using the SONICborden

This handpiece should be used in the same hose where the high-rotation handpieces (turbines) are used. When using it, the spray's air supply must be interrupted. For this, screw up either the screw (30), under the equipment's table depending on the hose installed where SONICborden is in use.

7.8.10 ULTRAscaler Use (accessory)

Remove the handpiece from the support and activate the key (8) of the foot control.

The frequency selector (61) is located on the table. In order to handle the ULTRAscaler, observe the instructions for use which come with it.

7.9 Operating the suction unit/ assistant's element

7.9.1 Handling of suction tubing

The suction system is actuated when the tubing is removed from the holder. To deactivate suction, replace tubing on holder.

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
7.9.2 Using the Poly 600 (accessory)

In order to handle the Poly 600 (24), observe the instructions for use which come with it.

7.9.3 Water jet for spittoon cleaning, through the register

The water jet to clean the spittoon is adjusted through the register (13) existing on the table.

7.9.4 Water jet for spittoon cleaning, with timed interruption, activated in the foot control (only for the Joystick foot control)

To obtain a flow of water to clean the basin, press the lever (11) or the key (56).

This flow will last for five, 15 or 30 seconds depending on the programmed time.

To regulate the water flow time, you must press and hold the lever for a while (11) or the key (56). The flow setup time is indicated by a beep.

- 1 Beep: 5 seconds of flow memorized;
- 2 beeps: 15 seconds of flow memorized;
- 3 beeps: 30 seconds of flow memorized;

Keep one of the functions listed above pressed until you hear the beep indicating the desired flow time.

The water volume is adjusted in the stopcock (13).

Exclusive for table 4T: The water supply is activated by activating the lever (15). It will remain activated between 5 and 60 seconds, according to the adjustment made in potentiometer (25). The water volume is adjusted in the stopcock (13).

7.9.5 Supply of the glass (accessory)

To get a water jet to supply the glass, press the driver (15) or (45).

7.9.6 Movement of the suction unit

The suction unit can turn 45°.

7.10 Operating the operating light

The operating light is an equipment with an intense, halogene light, appropriate for works which need lighting in the mouth area. Restorative resins are composites polymerisable by light and radiopaque. When using them we recommend to avoid the intense light of the operating light in the working field, so that it doesn't come to interfere with the polymerisation process.

We recommend using the appliance with the lamp positioned at a distance of approximately 70 cm from the operating field.

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Selection of the luminosity adjustment

2 adjustment modes of the luminous intensity can be programmed:

- 3 levels of luminous intensity
- continuous adjustment of the luminous intensity

The operating light leaves the factory programmed to work in the continuous mode.

In order to alternate between the modes:

- disconnect the appliance through key (9) or (46);
- press key (10) or (46) for at least 5 seconds;
- release key (10) or (46);
- press key (9) or (46);

a. Mode 1: 3 levels of luminous intensity

Viewing in the luminous indicator (38), on the table of the delivery unit or Pad Control (47)

1 LED on (or 1st) = minimum intensity $\approx$ 6 000 LUX

5 LEDs on (or 2nd) = medium intensity $\approx$ 13 000 LUX

All the LEDs on = maximum intensity $\approx$ 22 000 LUX

Turning on

Slightly press button (9) or (46) or with the chair at the patient's in-out-position, move it towards the working position.

The luminous intensity is the one adjusted upon the last disconnection of the appliance.

Alteration of the level of intensity

A sustained (long) pressure on key (9) or (46) alters the level of luminous intensity.

Example: if the appliance was disconnected at the intensity of 22 000 LUX.

A brief pressure turns the appliance on at the intensity of 22 000 LUX.

A sustained pressure alters the intensity to 13 000 LUX.

A new sustained pressure alters the intensity to 6 000 LUX.

A new sustained pressure alters the intensity to 6 000 LUX.

Turning off

Slightly press button (9) or (46) or move the chair to patient's in-out position.

During the chair's movements, the operating light will turn off, turning on again upon movement stop (except the patient's in-out position movement).

b. Mode 2: continuous adjustment of the luminous intensity

Viewing in the luminous indicator (38), on the table of the delivery unit:

10 levels of luminous intensity, from 1 LED on 1st LED on (6 000 LUX) to all the LEDs on (22 000 LUX).

Turning on

Slightly press button (9) or (46) or with the chair at the patient's in-out-position, move it towards the working position.

The luminous intensity is the one adjusted upon the last disconnection of the appliance.

Alteration of the level of intensity

Press key (9) or (46) and maintain it pressed. There will be a continuous alteration of the luminous intensity, in a descending order of value. Upon reaching the least value, the intensity returns to the maximum value. Another sustained pressure makes the intensity return to its maximum.

Turning off

Slightly press button (9) or (46) or move the chair to patient's in-out position.

During the chair's movements, the operating light will turn off, turning on again upon movement stop (except the patient's in-out position movement).

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



8 CLEANING / DISINFECTION

Upon utilization of disinfecting products, it is important to observe the usage prescriptions provided by the manufacturer of the products as well as the determinations of the public health agencies.

There are in the market several disinfecting products with different concentration levels. Upon their selection, take into account the following factors regarding the equipment:

- Nature of the surface to be disinfected and whether it may undergo corrosion or chemical attack.
- Contamination degree and its mode of elimination.
- Quality of the water used to dilute and its influence on the disinfection.
- Disinfection method.
- Safety upon handling and use.

We recommend disinfectants with a peracetic acid base (4 g of peracetic acid + 26 g of hydrogen peroxide for each 100 g of solution).

Disinfection shall be performed as follows:

- Wipe out the excess of contaminating matter by using either an absorbing paper or cloth;
- Discard this material into an adequate place;
- Apply the suitable disinfectant agent on the area to be treated, and let it there as longer as necessary;
- Wipe out the disinfectant agent with a wet cloth.

KaVo does not be held responsible for damages caused by:

- not recommended products/concentration levels
- products admittedly corrosive as sodium hypochlorite usually employed on floor tiling cleaning, which due to such application may have contact with the equipment.

When cleaning avoid liquids penetrating inside the equipment, as this may damage it.

8.1 Table with guidelines for routine preservation

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O= observation

C: 3F handpiece

T: Check air and water passage. Blocked passage at the nozzles can be cleaned with the help of a needle.

F: daily.

O: after placing the nozzle and before starting the treatment, actuate the air key to assure dry air during the treatment.

C: 3F handpiece.

T: clean the external side using 70% GL alcohol or other disinfecting agent without having an ammonium or acid basis or chlorine.

F: after each treatment.

O: Do never immerse the 3F handpiece in disinfection baths.

C: 3F handpiece

T: lubricate o-rings with liquid vaseline, without removing them.

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F: semanal.

C: 3F handpiece nozzle.
T: sterilize.
F: after each treatment.
O: nozzles can be sterilized in a vapor autoclave up to 135 °C.

C: spittoon bowl.
T: wash. Do not use acid or abrasive detergents nor detergents that produce foam.
F: several times a day..
O: to with draw it, pull it vertically with rotating movements.

C: spittoon bowl strainer.
T: wash. Do not use acid detergents or abrasives detergents nor detergents that produce foam.
F: several times a day..

C: suction hoses.
T: disinfect the suction hoses with a disinfectant for PVC hoses. Suction slowly with the instrument from 1 to 1,5 liters of disinfecting solution.
F: daily.
O: manufacturer's instructions for use must be obeyed rigorously. Excessive doses can be lead to a an early aging of the hoses.

C: Internal cleaning of tubings and pipes of the suction system.
T: wash out the inner pipes by suctioning water.
F: after each procedure that uses corrosive chemicals like sodium hypochlorite, in order to reduce the residual effects of such elements.

C: suction tubings connections.
T: remove connections from tubings and clean them with a brush and running water.
F: daily.

C: Suction tubings strainer.
T: wash with running water.
F: dily and after cirurgical tratments.
O: disassemble the parts in order to acess the strainer.

C: Protection and mirror of operating light.
T: clean.
F: weekly.
O: clean with a dry flannel. Avoid cleaning them with liquid agents.

C: Supply air net.
T: drain
F: daily.

C: Chair.
T: in order void marks at points where pressure is exerted, do not leave any object on the upholstered parts.

8.2 Technical safety controls

Responsibility: user or person in charge of preserving the product.

- Before using the equipment, assure yourself of the safety of its functions and of its perfect condition.
- Don't proceed with work in case there are parts damaged.
- Don't stick anything to the instruments' tubings.
- Do not blow spray or blowing air into open wounds during surgery since this may present a risk of air embolism or cutaneous emphysema.
- The maximum load of 135 kg for the upward movement must not be exceeded.
- Maximum tray holder loads: 500 kg
- Remove burrs and diamond tips from rotating handpieces, after completing treatment.
- Don't place the scaler in its holder without the tip wrench – danger of injuries and infection.
- After placing the nozzle and before starting the treatment, actuate the air key to assure dry air during the treatment.

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- Due o stagnation, the water and air ducts in treatment units are to be flushed or blown out when the unit is operated for the first time and after periods of standstill (after the weekend, after a holiday, vacation, etc.).
- Replace the water tanks of the spray and assepto sys, when the validity date indicated on them expires.
- Disconnect the unit main switch (supply electrical net), daily, before leaving the office.
- Close the water general register (supply water net), daily, before leaving the office.
- Close the air general register (supply air net), daily, before the interruption of equipment operation for longer time.

8.3 Annual inspections

They should be requested from authorized technical assistance network. The expenses arising from these services are not covered by the warranty.

- Check/redo the adjustment of the air pressure propelling the instruments.
- Verify the instrument spray.
- Clean the metallic screen in the water input connection, in the distribution box.
- Change the filtering element of the air filter, in the distribution box.
- Check/redo the adjustment of the valves regulating the static pressure of air input and pressurization of the tanks.
- Check/adjust the water outflow to the glass.
- Check/adjust the water outflow to the equipment syringe.
- Check/adjust the outflow of air and water to the suction system syringe.
- Change the o-rings of the triple syringe pointer.
- Check/adjust the outflow of air and water to the spray.
- Verify the stability of the equipment.
- Verify the correct movement of the equipment components.
- Verify the stability and perfect mechanical status of the articulated arms.
- Verify the braking systems.
- Verify the stop systems, limiting the chair movements.
- Verify the perfect connection of hoses in terminals.
- Handpiece tubings must be visually checked in regular intervals and replaced if damaged.

8.4 Assepto Sys system (accessory)

The Assepto system has the function to disinfect the inner water tubing of the instruments hoses (except syringe) by means of disinfectant solution action through the tubing. The reservoir (20), capacity of 0,50 liter, is located under the suction unit table.

To feed it:

- Unscrew the reservoir from its connection under the suction unit table;
- Fill the reservoir with disinfectant solution. We recommend disinfectants with a peracetic acid base (4 g of peracetic acid + 26 g of hydrogen peroxide for each 100 g of solution);
- Thread the reservoir (bottle) to its connection again only with the hands, enough to lean the bottleneck on the sealing rubber. After putting back the reservoir, the heard noise is not caused by air leakage, but by the pressurisation of the reservoir.
- The excessive tightening causes damage to the sealing set.
- After refilling and putting back the reservoir, drive an instrument with spray for some seconds to allow the air to go out the tubings.

How to operate

- Remove the tubing to be disinfected from the holder;
- Uncouple the handpiece;
- Actuate button (33) located under the dentist's element table, left side. The solution will go through the tube. If the tube is full, rest it on the holder while the disinfectant is acting;
- Follow the same procedure with the remaining tubings;
- After using the Assepto Sys, at any solution indicated, actuate the foot lever to get spray for at least 10 seconds to rinse the tubings.

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

The antimicrobial efficacy can only be confirmed by means of Official Methodologies, established by competent Organs. Thus, it is not feasible to execute tests in the place of use, as there is a risk of obtaining results which are not valid.

In order to ensure the Disinfection Process, we suggest that the confirmation of the antimicrobial efficacy of the product in use is provided by presenting Reports of Antimicrobial Tests executed by Official Laboratories and Methodologies.

Water must remain in the reservoir when the Assepto Sys system is not being used.

Do never drive the valve of the Assepto Sys (under the equipment's table) with instruments in the holders.

9 FAILURE ELIMINATION

Here are corrective measures that can be accomplished by the user. If not sufficient to settle the problem, please contact the KaVo Technical Service Net.

KaVo keeps available in its Technical Service Net all necessary information for the qualified technicians of the customer to repair the parts of the equipment considered repairable.

- P - Nothing is working
- C - Main switch OFF
 - Fuse rupture
- S - Turn the main switch on
 - Check system inlet fuse. If necessary, substitute it. Kavo reference for orders 0. 223.2783

System inlet fuse is together with the general switch, in the distribution box.

To change it:

- power down the equipment
- take off cover (67) of the fuse box
- change the fuse (66)
- replace the lid

- P - No handpieces spray
- C - Spray not selected
 - Water tubing register, closed
 - Water main register, closed
 - Compressor off
- S - Select spray
 - Open tubing spray register
 - Open water main register
 - Turn compressor on

- P - Leak at tubing; handpiece / coupling
- C - Defective rubber rings at coupling
- S - Substitute O-rings

- P - Suction tubings do not operate properly
- C - Sliding valve at tubing connection closed
 - Blocked filters
 - Blocked suction tubings
 - Suction pump does not work
- S - Open sliding valve
 - Substitute filters
 - Wash tubings with a lot of water
 - Check suction pump

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10 OPTIONS AND ADAPTABLE ACCESSORIES

Options and accessories listed below apply exclusively.

Options

- Table T:
A table with hanging hoses

- Table S:
A table with hoses adapted in moveable arms

- Table T5
A table with hanging hoses5

Accessories

- Children's seat

Upholstered seat and backrest suited upon the normal upholstery, for use in Odontopediatrics.

- Fourth instrument

This kit provides all necessary mechanical and hydropneumatical parts to allow the use of a fourth handpiece.

- Fourth instrument multiflex LUX

This kit provides all necessary electrical, mechanical and hydropneumatical parts to allow the use of a fourth handpiece Multiflex LUX.

- Installation kit for vacuum pump

This kit provides all necessary mechanical, hydropneumatical and electrical parts to allow the use of a vacuum pump. The pump is not part of the kit.

- Assepto Sys

This kit provides the asepsis system of the handpieces tubings through passage of a disinfectant solution.

- Water cock for the glass

This kit provides all mechanical and hydraulics parts for the installation of a water cock for the glass.

- Kit Air inlet

Air filter with automatic drain to be fitted to the equipment's junction box

- Autonomous distribution box kit

Autonomous cover interconnected to the box integrated with the chair, in order to allow the installation of the unit at places where the connection spot is already prepared, and does not comply with the standard positioning.

- Kit water filter

This kit has the function to supply water to the tap on suction units

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

- Standard Headrest

Double-jointed Headrest

- Kit Auxiliary arm for the suction unit

Composed of two jointed arms and one handpiece holder and applied to the column of the suction unit. Replaces the suction unit's handpiece holder, which is deactivated.

- ULTRA scaler

Ultrasonic appliance for application in Prophylaxis.

- Tray with a rim:

A stainless steel tray to hold instruments (used on the dentists table).

- Auxiliary Syringe:

The kit consists of the mechanical and hydropneumatic equipment necessary to use a triple syringe.

- Venturi Suction Kit-Ø 9.5:

The kit consists of the mechanical and hydropneumatic equipment necessary for operating the suction hose for 9.5 mm. diameter dental suction cannula's.

- Borden tubing Set

Universal Tubing with Borden coupling.

- Midwest tubing Set

Universal Tubing with Midwest coupling.

- Multiflex tubing Set

Universal Tubing with Multiflex coupling.

- EXHAUST Set

(Only Midwest tubing)

- Keyboard to S table

(Keyboard with functions + left pad without functions)

-Keyboard to Scaler

(Left Keyboard with functions + protections)

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- Table + Adjustable arm

The table is attached to rotating a fixed to the suction nozzle support. The chair follows the height when it is in movement. The height of the column can be adjusted.

- Table + arm (pantograph with pneumatic brake)

The table is attached to rotating arms with jointed, fixed to the suction nozzle support. Follows the height of the chair when in movement. Some predetermined height can be adjusted manually.

- Table + arm (pantograph)

The table is attached to rotating arms with jointed, fixed to the suction nozzle support. Follows the height of the chair when in movement. Some predetermined height can be adjusted manually.

- Table and movable stand (Cart)

The table fits a pedestal. It has four rollers to allow extra mobility. The height of the table is fixed.

- Ortho table

The table fits a rotating arm system fixed to the column of the suction device. It was developed to help orthodontist activities.

- Comfort upholstery**- Standard upholstery****- 5th holder and 6th holder Kits**

Holders for 3F handpiece or ULTRAscaler, fixation in the dentist's table.

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

Suportes a serem fixados na mesa standard, para seringa triplíce ou ULTRAscaler

11 SIMBOLOGY

	Applied part type B
	Applied part type BF
	Caution, refer to accompanying documentation
	Chair down movement
	Chair up movement
	Backrest forward movement
	Backrest rearward movement
	Position in-out-of patient movement
	Last Automatic Position
	In-out of Patient Position
	Automatic movement for a first work position
	Automatic movement for a second work position
	Backrest forward and rearward movement
	Fuse
	Main switch - operating light
	Polymerizator
	Alternating current
	Water for spittoon bowl cleaning
	Tumbler water
	Cooling spray
	Assepto Sys
	Foot control of the dentist's element
	Handpiece with light
	Turbine
	3F handpiece
	Continuous function with intermittent load
	Electrical MicroMotor
	On
	Off

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Electric supply cable: 3 x 2.5mm² / 500mm

Input voltage/frequency:

The equipment is default adjusted for operation with 220V~net voltage.

However, it may be changed to operate with net voltage = 110V~ or 127~ / 50/60 Hz

Max. power consumption: 280 VA

Inlet fuse:

110V~: T 4 A

127V~: T 2,5 A

220V~: T 1,6 A

Water connection: R1/2"

Water inlet pressure: 2.5 – 6 bar

Water hardness in millimole by liter:

1) Total hardness: 1,074 – 1,790

2) Carbonate hardness: 0,716 – 1,074

Drain connection in mm: Ø 40

Pipes flow: 5 l/min.

Air connection: R1/2"

Air inlet pressure: 5 – 5,5 bar

Air inlet flow: min 100 l/min

Suction connection in mm: Ø 40

Spray water reservoir capacity: 1000 ml

Assepto Sys reservoir capacity: 500 ml

Luminous intensity of the operating light: 6 000 a 22 000 LUX

Unit – net weight (without accessories): T: 134 Kg S: 136 Kg

Maximum acceptable load of additional equipment on the chair: 60 kg

Maximum capacity of raising the chair (with max. load of additional coupled equipment): 135 kg

Maximum torque capacity of the accessory connection points: additional weights are not foreseen for this product.

Load capacity of the double-articulated headrest: 10 kg

No other additional weights are foreseen for the other accessories.

Colour codes of the conducting tubes:

The tubes which conduct air have blue markers and those conducting water are marked green.

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13 INSTALLATION

-Important notes

The installation services must be executed only an authorized technician, duly identified.


-Prerequisites for proper operation

The safety, reliability and performance of the equipment are warranted subject to the following conditions:

- The electrical installation of the surgery must comply with the requirements and specifications.
- The unit must be operated in accordance with the operating instructions.
- Assembly, extensions, adjustments, modifications or repairs must be carried out by technicians authorized by KaVo, and the completion of the installation must be reported to KaVo by means of a delivery record

Transport and storage

Follow the instructions printed on the outside of the packaging which are applicable to transport and storage.

Preparation of the place to install the product

The preparation of the installation area by the customer beforehand should follow the "Instructions for preparation of the installation place, for KaVo Unik".

This document contains information on the positioning and form of tubing and supply values.

- Covers for above floor connection are not included in the delivery.

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Terms of Warranty

§1. The product described is guaranteed against manufacturing defects within a period of 12 months from the effective date of installation of the product, except for upholstered parts for which the warranty period is 6 (six) months.

The warranty period starts from the date of installation by KaVo authorized technician provided that the product has not been stored for more than 03 months from the date of the issue of invoice (the eventual storage of the product will be to the account and risk of the customer).

§2. In cases where the storage period exceeds the period highlighted in the previous paragraph, the warranty period will start from the date on which that period expired.

§3. In case of a special extension of the warranty period, the latter does not cover the rubber made components such as, for instance, rings, diaphragms, sealings, and tubing.

§4. Lamps and fuses isn't included in the warranty.

§5. The warrantee is limited to the repair or replacement of parts with a manufacturing defect, duly noted by the manufacturer or Authorized Service Network, and does not including the repair of defects caused by:

- Disregard of the use and maintenance instructions
- Falls or bumps
- Inadequate storage
- Action of natural agents
- Installation and technical service carried out by personnel that are not authorized by KaVo
- Damages to the upholstered parts caused by clothes (stains), cutting or sharp objects and dental material
- Damage to parts painted and painted plastic parts or not, caused by misuse, dental materials such as acrylic, gloves, etc., and sharp or pointed objects
- Connecting the unit to the power network having a voltage different to that specified for the equipments operation
- Lack of grounding or inefficient grounding
- Use of accessories, spare parts and materials not authorized by KaVo.

§6. The warranty will be automatically invalidated:

- For the normal course of the period of validity
- The introduction, by the customer, of changes to the equipment or improper use of accessories
- Adulterations to the warranty certificate.
- For installation and technical assistance made by persons that are not authorized by KaVo
- For the storage for more than 15 months from the date of issue of the invoice

§7. Repair or replacement of parts during the warranty period does not extend the period of validity of the original warrantee.

§8. Parts replaced due to the guarantee will be the property of KaVo do Brasil Ind. e Com Ltda.

§9. The buyer will have to pay all duties that are due (regardless of the product warrantee):

- Installation of product (which will be informed at the time of the sale of the product);
- The cost of the visit for technical support from the seventh month of the warranty period;
- The travel and subsistence of staff involved in answering the call for the installation or technical support, based on the value of distance traveled (km) to be informed in the Customer Support Network KaVo.
- The labor, materials, parts and necessary adequations for site preparation for the installation of the product, mains power, air, water, drainage, suction, grounding and etc.

§10. KaVo do Brasil Ind. e Com Ltda. does not authorize any person or entity to carry out on its behalf, any work relating to the guarantee of their products than those specified here.

§11. The buyer must sign the delivery protocol of the product for installation. At the end of each visit he must also sign the technical support report with the services performed, parts replaced, etc., registering his approval. The refusal of the buyer to sign this report does not form an allegation of non-compliance.

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